

PHARMACEUTICAL CAPSULES

ACGcaps™ GF

G GELATIN



Make it better.

Our fish gelatin capsule.

ACGcaps™ GF derives from fish by-products and therefore contributes to environmental sustainability. These capsules are also a safe and globally accepted alternative to mammalian-sourced gelatin. Due to their low WVTR and OTR, they are suitable for moisture-sensitive and oxygen-sensitive drugs.

Key features

- / Low WVTR (water vapour transmission rate) and OTR (oxygen transmission rate)
- / TSE / BSE free
- / Suitable for moisture-sensitive and oxygen-sensitive drugs
- / Proven machine-friendly design



Technical specifications

- / Moisture content: 13.0% to 16.0%
- / Disintegration time in water at $37 \pm 2^{\circ}\text{C}$: maximum 15 minutes (with guided disc)
- / WVTR: $819 \text{ gm} / (\text{m}^2 \cdot \text{Day})^*$
- / OTR: $6.57 \text{ cc} / (\text{m}^2 \cdot \text{Day})^*$

* When exposed to 40°C / 75% RH